

Superior Testing at an Affordable Price

3300 Plastics Testing Package

The Instron® Model 3300 Plastics Packages are complete, out-of-the-box mechanical testing instruments for tension and flexure testing of plastic materials. With these packaged systems, you can conduct quality control or component tests to ensure your manufacturing processes meet the high standards your customers expect.

Key Features

- Method-driven for simple or sophisticated tests
- A large selection of additional grips and fixtures are available to meet future testing needs
- Extensometers available for strain measurements
- World-class precision and accuracy
- Excellent long-term investment
- Safe and reliable
- Quick delivery

Bluehill® Software

Operators can setup a test using one of the pre-configured test method templates, gather results using the extensive calculations library, and generate an email report when the test is complete.

To begin a test simply open a method and press start – it's that easy. With multiple operators, varying skill levels, and continuous training needs, you want a user interface that is simple to learn. Bluehill provides all the capabilities you need to handle basic test requirements quickly and efficiently. Test control, data acquisition, plotting, calculations, and reporting are performed automatically. To further assist the operator, the Bluehill prompted method feature guides operators through repetitive test procedures with step-by-step instructions for all user levels.



Plastics Testing Packages include Bluehill Software, PC, Installation, and Training. Systems can be **configured** with a variety of grips, flex fixtures, and extensometers **to meet ISO 527, ASTM D790, ASTM D638, and more.**



For information on this and other packages available, visit go.instron.com/3300packages or call us at 1.800.564.8378